


## Bos d 8

Native allergen protein

Catalog # **BD08-A9010**

Lot # 1Y2253-5

<b>Catalog Number</b>	<b>BD08-A9010</b>
<b>Allergen</b>	<b>Bos d 8</b>
<b>Alternate Names</b>	Alpha-S1-casein, Alpha-S2-casein, Beta-casein, Kappa-casein
<b>Source</b>	<i>Bos taurus</i>
<b>Accession Numbers</b>	Genbank: <a href="#">NP_851372</a> , <a href="#">NP_776953</a> , <a href="#">XP_005902099</a> , <a href="#">NP_776719</a> UniProt: <a href="#">P02662</a> , <a href="#">P02663</a> , <a href="#">P02666</a> , <a href="#">P02668</a>
<b>Host</b>	Native protein
<b>Formulation</b>	Native protein stored in 1X PBS solution.
<b>Molecular Weight</b>	Calculated: 20-30 kDa Observed: ~30 kDa
<b>Purity</b>	The purity of Bos d 8 was determined to be <b>&gt;80%</b> by densitometry.
<b>Concentration</b>	5.64 µg/µl



### Product Description

Native *Bos taurus* allergen.

### Storage and Stability

Store at -20°C. To avoid repeated handling and multiple freeze/thaw cycles aliquot diluted product into smaller quantities.

### References

1. Stewart AF, et al: Complete nucleotide sequences of bovine alpha S2- and beta-casein cDNAs: comparisons with related sequences in other species. *Mol Biol Evol.* May;4(3):231-41, 1987.
2. Stewart AF, et al: Nucleotide sequences of bovine alpha S1- and kappa-casein cDNAs. *Nucleic Acids Res.* May 11;12(9):3895-907, 1984.

To place your order, please contact us by phone 1-778-326-0223 or 1-888-606-3424 (Toll free) or by email: [orders@signalchemdx.com](mailto:orders@signalchemdx.com) or [info@signalchemdx.com](mailto:info@signalchemdx.com) - [www.signalchemdx.com](http://www.signalchemdx.com)

**FOR IN VITRO RESEARCH PURPOSES ONLY. NOT INTENDED FOR USE IN HUMANS OR ANIMALS.**

# SAFETY DATA SHEET

## Article 1 – Product Identification

**Product Name: Bos d 8****Catalog # BD08-A9010**

*This product is sold only for research use by qualified laboratory personnel, and is not to be used as a drug, medical device, food additive, cosmetic, nor household chemical. It is not to be used in diagnostic, therapeutic, consumer, agricultural, nor pesticidal applications.*

Supplier of Datasheet: SignalChem Diagnostics Inc.  
Street Address: 190-13160 Vanier Place  
City, Prov. Postal Code: Richmond, BC, V6V 2J2  
Country: Canada  
Emergency Phone: 1-888-606-3424 (Toll free)  
1-778-326-0223 (local)

## Article 2 - Hazard Identification

- **WHMIS Classification:** Not WHMIS controlled.
- **GHS classification:** Not GHS classified.
- **Hazard Pictograms:** No labelling applicable.
- **Signal words:** None.
- **Hazard statements:** None.
- **Precautionary statements:** None.
- **Other hazards:** None known.

## Article 3 – Composition/Information on Ingredients

**Chemical Characterization:** Mixtures.

**Description:** This product consists of the substances listed below.

Common name	Chemical name	CAS-No.	Concentration
Sodium Phosphate, Dibasic	Sodium Phosphate, Dibasic	7782-85-6	0.248%
Potassium Phosphate, Monobasic	Potassium Phosphate, Monobasic	7778-77-0	0.024%
Sodium Chloride	Sodium Chloride	7647-14-5	0.72%
Potassium Chloride	Potassium Chloride	7447-40-7	0.02%
Protein	N/A	N/A	≤0.1%

## Article 4 – First-aid Measures

- **General information:** Consult a physician by providing the SDS.
- **After inhalation:** In case of irritation by inhaling this product, move affected person to fresh air and await recovery. If irritation persists, seek immediate medical attention. If casualty cannot breathe, give artificial respiration and seek immediate medical attention.
- **After skin contact:** Immediately wash with soap and plenty of water and rinse thoroughly. Consult a physician.
- **After eye contact:** Rinse opened eyes with plenty of water for at least 15 minutes. Consult a physician.
- **After swallowing:** Rinse the mouth with plenty of water and consult a physician.

## Article 5 - Fire-fighting Measures

- **Suitable extinguishing media:** Use water spray, extinguishing powder, carbon dioxide, or other appropriate measure that is suitable to the environment.
- **Specific hazards arising from the substance or mixture:** None known.
- **Special protective equipment and precautions for fire-fighters:** Self-contained breathing apparatus if necessary.

## Article 6 – Accidental Release Measures

- **Personal precautions, protective equipment, and emergency procedures:** Apply standard laboratory practices and personal protective equipment. Avoid breathing vapors, mist, or gas. Ensure adequate ventilation.
- **Environmental precautions:** Do not allow to enter drains.
- **Methods and materials for containment and cleaning up:** Absorb on sand or vermiculite and place in closed containers for disposal.

To place your order, please contact us by phone 1-778-326-0223 or 1-888-606-3424 (Toll free) or by email: [orders@signalchemdx.com](mailto:orders@signalchemdx.com) or [info@signalchemdx.com](mailto:info@signalchemdx.com) - [www.signalchemdx.com](http://www.signalchemdx.com)

**FOR IN VITRO RESEARCH PURPOSES ONLY. NOT INTENDED FOR USE IN HUMANS OR ANIMALS.**

# SAFETY DATA SHEET

## Article 7 - Handling and Storage

- **Precautions for safe handling:** Wear chemical safety goggles and compatible chemical-resistant gloves. Avoid inhalation, contact with eyes, skin, or clothing.
- **Conditions for safe storage:** Store according to product label instructions. Keep container upright and tightly closed.

## Article 8 - Exposure Controls/Personal Protection

- **Components with limit monitoring values at workplace:** N/A
- **Appropriate engineering controls:**  
Apply adequate ventilation including mechanical exhaust or laboratory fume hood. Follow standard laboratory practices.
- **Individual protection measures:**
- **Respiratory protection:**  
Use appropriate respirator if there is inadequate ventilation by following the government standards.
- **Hand protection:**  
Wear gloves and use proper glove removal technique to avoid skin contact. Discard gloves after use by following the applicable laboratory regulations. Wash and dry hands.
- **Eye/face protection:**  
Safety goggles with side-shields approved under appropriate government standards.
- **Skin/body protection:**  
Use appropriate clothing, footwear and any additional protection measures to protect from splashing or contamination.

## Article 9 – Physical and Chemical Properties

<b>Appearance:</b> Colorless liquid.	<b>Danger of explosion:</b> Not determined.
<b>Odour/Odour Threshold:</b> Not determined.	<b>Explosion limits:</b> Not determined.
<b>pH:</b> Not available.	<b>Decomposition temperature:</b> Not available.
<b>Melting point/freezing point:</b> Not determined.	<b>Vapor pressure at 20 °C:</b> Not determined.
<b>Boiling point/Boiling range:</b> Not determined.	<b>Density:</b> Not determined.
<b>Flash point:</b> Not determined.	<b>Relative density:</b> Not determined.
<b>Flammability (solid, gaseous):</b> Not determined.	<b>Vapor density:</b> Not determined.
<b>Ignition temperature:</b> Not determined.	<b>Evaporation rate:</b> Not determined.
<b>Auto-igniting:</b> Not determined.	<b>Solubility in / Miscibility with Water:</b> Fully miscible.

## Article 10 - Stability and Reactivity

- **Reactivity:** Stable under recommended transport and storage conditions.
- **Chemical stability:** Stable under recommended transport and storage conditions.
- **Possible hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:** Heat and moisture.
- **Incompatible materials:** Strong acids/bases, strong oxidizing/reducing agents.
- **Hazardous decomposition products:** Carbon oxides may be formed under fire conditions; no known decomposition information for other decomposition products.

To place your order, please contact us by phone 1-778-326-0223 or 1-888-606-3424 (Toll free) or by email: [orders@signalchemdx.com](mailto:orders@signalchemdx.com) or [info@signalchemdx.com](mailto:info@signalchemdx.com) - [www.signalchemdx.com](http://www.signalchemdx.com)

**FOR IN VITRO RESEARCH PURPOSES ONLY. NOT INTENDED FOR USE IN HUMANS OR ANIMALS.**

# SAFETY DATA SHEET

---

## Article 11 - Toxicological Information

---

- **Acute toxicity:** Not available.
- **LD/LC50:** Not available.
- **Skin corrosion/irritation:** Not available.
- **Serious eye damage/eye irritation:** Not available.
- **Respiratory or skin sensitization:** Not available.
- **Germ cell mutagenicity:** Not available.
- **Carcinogenicity:** No components are listed in IARC, or NTP, or OSHA, or ACGIH.
- **Reproductive toxicity:** Not available.
- **Teratogenicity:** Not available.
- **Specific target organ toxicity - single exposure/ - repeated exposure (GHS):** Not available.
- **Aspiration hazard:** Not available.
- **Potential health effects:**
  - Inhalation:** Not available.
  - Ingestion:** Not available.
  - Skin:** Not available.
  - Eyes:** Not available.
- **Signs and Symptoms of Exposure:** Not available.
- **Synergistic effects:** Not available.

---

## Article 12 - Ecological Information

---

- **Eco-toxicity:** No data available.
- **Biodegradability:** Not applicable.
- **Bio-accumulative potential:** Not applicable.
- **Mobility in soil:** Not applicable.
- **PBT and vPvB assessment:** Not applicable.
- **Other adverse effects:** Not applicable.

---

## Article 13 - Disposal Considerations

---

- **Disposal methods:** In accordance to applicable national, regional, or local laws and regulations. For additional handling information and protection of employees please refer to Article 7 and 8.
- **Contaminated packaging:** Disposal should be made in accordance to official regulations. Use water or cleansing agents to clean the area.

---

## Article 14 - Transport Information

---

- **DOT:** Not dangerous goods.
- **IMDG:** Not dangerous goods.
- **IATA:** Not dangerous goods.

---

## Article 15 - Regulatory Information

---

- **WHMIS Classification:** Non-hazardous.
- **GHS label elements:** Not applicable.
- **Signal word:** Not applicable.
- **Hazard statements:** Not applicable.

---

## Article 16 - Other Information

---

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. SignalChem shall not be held liable for any damage resulting from handling or from contact with the above product. See the Technical Specification, Packing Slip, Invoice, and Product Catalog for additional terms and conditions of sale.

To place your order, please contact us by phone 1-778-326-0223 or 1-888-606-3424 (Toll free)  
or by email: [orders@signalchemdx.com](mailto:orders@signalchemdx.com) or [info@signalchemdx.com](mailto:info@signalchemdx.com) - [www.signalchemdx.com](http://www.signalchemdx.com)

**FOR IN VITRO RESEARCH PURPOSES ONLY. NOT INTENDED FOR USE IN HUMANS OR ANIMALS.**